

PSME2-proteasome (prosome, macropain) activator subunit 2 (PA28 beta) Gene

Reference: PTXBC072025

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PA28B; PA28beta; REGbeta; proteasome activator complex subunit 2; 11S regulator complex beta subunit; 11S regulator complex subunit beta; MCP activator, 31-kD subunit; REG-beta; activator of multicatalytic protease subunit 2; cell migration-inducing prote		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) activator subunit 2 (PA28 beta)
GENE ID	5721
ORIGIN SPECIES	Human
ACCESSIONS	BC072025
VECTOR	pDONR223
NCBI SYMBOL	PSME2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaagccgtgtgggtgcgcctgagcggggaagccgcaaacagggtggaggtcttcaggcagaatctttccaggagcctgaggaaatcctctacagattctgccacagaaaatcatatacctgaatcagctcttgcagaggactccctcaatgtggctgacttgacttcctccg
ggccccactggacatccccatccagacccctccaccaaggatgatgagatggaaacagataagcaggagagaagaagagtcctaaagtgtggatttctccctgggaatgagaaagtcctgtccctgcttgcctgtgtaagccagaagctggactctcaagagaaatgcattctggtg
attacatggatccaacacctgatcccaagattgaagatggaaatgattttgggtagcaatccaggagaaggtgctggagaggtggaatgccgtcaagaccaaagtggagccttccagacaacattccaagtactctcagaacgtgggatgctgtggccaaggcctccaaggaga
ctcatgtaatggattaccgggcttggtgcatgagcgagatgaggcagcctatggggagctcagggccatggtgctggacctgagggccttctatgctgagctttatcatatcatcagcagaacctggagaaaattgtaacccaagggtgaagaaaagccatctatgtactga