

OXCT2-3-oxoacid CoA transferase 2 Gene

Reference: PTXBC072388

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FKSG25; SCOTT; succinyl-CoA:3-ketoacid coenzyme A transferase 2, mitochondrial; 3-oxoacid CoA-transferase 2A; testis-specific succinyl CoA:3-oxoacid CoA-transferase; 3-oxoacid CoA-transferase 2		

For research use only.

GENE DESCRIPTION	3-oxoacid CoA transferase 2
GENE ID	64064
ORIGIN SPECIES	Human
ACCESSIONS	BC072388
VECTOR	pDONR223
NCBI SYMBOL	OXCT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgctgcggctcctggtcgtcagtcggcgccgggtccccggcggtcagggctcgcgtgtccagggtgcgccgctgtttgccaccagtccccggccccgtccaagtctacgcggaccgggtggagatggtgaaggacatctctgacggcgaccgt
catgatcggggcttcgggctctgctgggatccccgagaacctgatcgcgcgtgtcagggaccgcgtgaaagacctgcaggtgtcagcagcaacgtggcgctggaggaattcggcctcctgctggcccgccaggcaggtccgtcgcacgtctgttctacgtggcgga
gaacacctgtgcgagagccagttacgtggcaggagagctggagctggagctcagccccaggcgaccctggccgagcgcacccgcggggcgccgggtgcccgcctctacacccccagggtacgggaccctggtccaggagggggcgccccatccgctacacc
ccggacggccacctggcgctcatgagccagccccgagaggtgaggagttcaacggcgaccacttcttttggagcgccatccgggcagacttgcctgtgtgaaagggtggaaggccgaccggcgaggaaacgtgtgtctcaggagagcgccgcaatttcaacgtgccat
gtgcaagctgcagacgtcacggcggtggaggtgggggtctccccccagaagacatccagttcctaacttatgtag