

C11orf57-chromosome 11 open reading frame 57 Gene

Reference: PTXBC072455

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C11orf57; chromosome 11 open reading frame 57		
For research use only.			

GENE DESCRIPTION	chromosome 11 open reading frame 57
GENE ID	55216
ORIGIN SPECIES	Human
ACCESSIONS	BC072455
VECTOR	pDONR223
NCBI SYMBOL	C11orf57
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgaaaataagagaactggaataacagatggaagatgctaccgggggacaaaaggaaaatgctaccagcagttcaagccgatgcgcagtgatggttttgatgaagaaagtcaagatactattggaggccaaagaatgaaattctgggacactggaagatgattttcttaagg
ctaaatcctggaataaaagtctatgattatgaagcaaacatgccagacagatgggtcacagtggttataaagagtataacctgaagaatttgaacagacagtgatcagcaagatattaccaacgggaagaaaacatctccccaggtaagtcctaccatgaatccgcaaacaca
agaagtcaagaaatcccaaaaaaagcagaaaaaagggtcacacaaaaacagaagaaagcaaaaaggaaagccacagatataacagcagatcctcagtgagttctcagaagaaactggggtcttggtacaaaggaaaggaaacacacataaacgcaagaaaaatcc
aggaaaaagtctcctcaaaaacctgctttattcttagggcagaagtaacacttcacattcagatgattcagcatccagcagttctgaggaaagtgggaaagagacactaagaaacaaaaggaaaaagagagaaaaagcccatcctctgtgccaacaatgaatacaggag
aggacaaacaaacgcacaaattggaaagtagctacagatgaaggtctgctgagagctcagaggatgactaa