

UBAC2-UBA domain containing 2 Gene

Reference: PTXBC072674

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PHGDHL1; ubiquitin-associated domain-containing protein 2; RP11-178C10.1; UBA domain-containing protein 2; phosphoglycerate dehydrogenase like 1; phosphoglycerate dehydrogenase-like protein 1; UBA domain containing 2		

For research use only.

GENE DESCRIPTION	UBA domain containing 2
GENE ID	337867
ORIGIN SPECIES	Human
ACCESSIONS	BC072674
VECTOR	pDONR223
NCBI SYMBOL	UBAC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccggtctgtgtctacagacagcaaaatgtccaggtgcatcaggtgctctgcatccccagctggatggcaaaattcttttggacacttgaacccatcttctcttcagaacccaccagcgaagccagaattgggatgggagccacgtggacatccagagacagagagaatggag
ctgtggaccggcagctgatgttctcagtttgcacaaggaggcgacagagacagcagcagggagggaatgatcaattggaatcgtcttttctctctttacgtcagcagacaaaacgtaaactatcaggcggtcgcagctctgagccagcagcgcccccctagaagtctctgaggaac
aggctgccggctcatggagatgggattttccagaggtgatgctttggaagccctgagagcttcaacaatgacatcaatgtcgccaccaactctctgctgagcactga