

C2orf79-chromosome 2 open reading frame 79 Gene

Reference: PTXBC073803

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C2orf79; peptidyl-tRNA hydrolase domain-containing protein 1; peptidyl-tRNA hydrolase domain containing 1		

For research use only.

GENE DESCRIPTION	chromosome 2 open reading frame 79
GENE ID	391356
ORIGIN SPECIES	Human
ACCESSIONS	BC073803
VECTOR	pDONR223
NCBI SYMBOL	C2orf79
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacccggggagtaggtccggccttcgggtgtcaggaagatggcgcctctggggcggagccgaggtctgtgacatactgtgtgtacgaaaggatctatcacaaagctcgttctctgcccggcgggcgactggtagcgcaggtgtcacgcggccaccgcggccttgc
acaactcacccgcaccacccgcacacagccgcttacctccaagagctggggcgcatgcgcaaaagtgtcctcgaggccccagatgagaccacctaaaggagctgcccagaccctgcaacagaagaacattgaccacatgctgtgcttgagcaaccagagaatatcgccacttga
ttgctctccggccctacccaaggagaagtgggccagtatttgaagaagttccgattgttcaagtaa