



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZFC3H1-zinc finger, C3H1-type containing Gene

Reference: PTXBC073843

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

CCDC131; CSRC2; PSRC2; zinc finger C3H1 domain-containing protein; coiled-coil domain containing 131; coiled-coil domain-containing protein 131; proline/serine-rich coiled-coil 2; proline/serine-rich coiled-coil protein 2; zinc finger C3H1-type containing

For research use only.

GENE DESCRIPTION	zinc finger, C3H1-type containing
GENE ID	196441
ORIGIN SPECIES	Human
ACCESSIONS	BC073843
VECTOR	pDONR223
NCBI SYMBOL	ZFC3H1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr