

FLJ32065-hypothetical protein FLJ32065 Gene

Reference: PTXBC073870

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein FLJ32065
GENE ID	201283
ORIGIN SPECIES	Human
ACCESSIONS	BC073870
VECTOR	pDONR223
NCBI SYMBOL	FLJ32065
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaaatagtagcggtattccgaacagacactaaaaatagctgtcatctcaagaatccagtgcctgtgtcacagtatgagaaagtagatgtggggaacagcgtttaatgaatgaagcattccagccagcagtgatctcttggacctgcattctccatcagattggatcacctcccacctg
aggcccccaagactttgaacagttcttcagtcattcattacagaagaatccctctccagacaaacgcagtaattatatacgggtccattggatctctatgaagcaccagaattatcagtgagaataatataaaggctcacgggctactgtaaacatattctatgccttga