



Reference: PTXBC075814

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRIPR; PGI receptor; PGI2 receptor; prostaglandin I2 receptor; prostanoid IP receptor; prostaglandin I2 (prostaglandin) receptor (IP)		

For research use only.

GENE DESCRIPTION	prostaglandin I2 (prostacyclin) receptor (IP)
GENE ID	5739
ORIGIN SPECIES	Human
ACCESSIONS	BC075814
VECTOR	pDONR223
NCBI SYMBOL	PTGIR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]