

CC2D2B-coiled-coil and C2 domain containing 2B Gene

Reference: PTXBC075861

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein CC2D2B; C10orf130; bA248J23.4; coiled-coil and C2 domain containing 2B		
For research use only.			

GENE DESCRIPTION	coiled-coil and C2 domain containing 2B
GENE ID	387707
ORIGIN SPECIES	Human
ACCESSIONS	BC075861
VECTOR	pDONR223
NCBI SYMBOL	CC2D2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttcccaaccgaagaatcgtaactactgttttaataatgaaggatagcttctagtcacagatataatcaaggcattaaatccacctcagcaactctggatataatttccacaattctaatgcaacatttgacctcattgtctgattgtatctttgattcctttgtgcctaatacaccagatgaa
aatgatggctctgatatatggatgacatcagagcactgcatcagtttagctatcggaataagggagagcatgccatccttctgtgaattcttctgtattttggaagaaggcactggtcctctgggaacgtcagtgtagaaggcatgtggcttatgtagtaactcaagaactaatgaatat
ttgctttggaatcacaactggccaatgttataagcagtttgacccgtttgtcccttaaaagtgtagattgtttgttgatgatagaatgtctgtgttaataattcaacaaaataatacaccatggctgtatttttactattcaaggaaagtgtctggaagcagttgctccaaaaacgttcaagg
aacaaaaatacaaaagcatacaggtcagggatftcccatccagatgccatacattgatgtacagtcattattgatgtctgtttatcaaaactggaattcactctgtgaatttccccagacagaattgtcttagctgtatactaccatacccaacaacataattatctgtgtgggtctatttggttc
cttagtcaacatcaatga