



Reference: PTXBC077718

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CK13; K13; WSN2; keratin, type I cytoskeletal 13; CK-13; cytokeratin 13; keratin 13, type I; keratin 13		

For research use only.

GENE DESCRIPTION	keratin 13
GENE ID	3860
ORIGIN SPECIES	Human
ACCESSIONS	BC077718
VECTOR	pDONR223
NCBI SYMBOL	KRT13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]