

DAAM2-dishevelled associated activator of morphogenesis 2 Gene

Reference: PTXBC078153

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dJ90A20A.1; disheveled-associated activator of morphogenesis 2; dishevelled associated activator of morphogenesis 2		

For research use only.

GENE DESCRIPTION	dishevelled associated activator of morphogenesis 2
GENE ID	23500
ORIGIN SPECIES	Human
ACCESSIONS	BC078153
VECTOR	pDONR223
NCBI SYMBOL	DAAM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccccgcaagaggagccaccatggcctgggtctctgtgctgttcggggcagtgacatccccgaaatcaacctccgggacaaccacctctgcagttcatggagtttccagccccatctcgaacgcagaggagctcaacatccgcttgcagagctggatgaattggat
ctcactgacaaaaccgagaggctatgtttgcactgccccctgagaaagaatggcagatctactgcagcaagaagaagggtgccctctctgacccctctgcccacatctcagggttcagtgcatgggtggctcttgccttggcctgctcatgcattcacctgatgtttattacatgccagc
cctgttctaggtgttgagggaacaacagtgaataa