



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM35-transmembrane protein 35 Gene

Reference: PTXBC078658

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

TMEM35; transmembrane protein 35A; transmembrane protein 35

For research use only.

GENE DESCRIPTION	transmembrane protein 35
GENE ID	59353
ORIGIN SPECIES	Human
ACCESSIONS	BC078658
VECTOR	pDONR223
NCBI SYMBOL	TMEM35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]