

GMFG-glia maturation factor, gamma Gene

Reference: PTXBC080180

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GMF-GAMMA; glia maturation factor gamma		

For research use only.

GENE DESCRIPTION	glia maturation factor, gamma
GENE ID	9535
ORIGIN SPECIES	Human
ACCESSIONS	BC080180
VECTOR	pDONR223
NCBI SYMBOL	GMFG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgactccctgggtggtgtagagtagacccagagctaacagaaaagctgaggaaattccgttccgaaaagagacagacaatgcagccatcataatgaagggtggacaaagaccggcagatgggtgctggaggagaatttcagaacatttcccagaggagctcaaatgga
gttgcggagagagacagcccaggttcgtgtttacagctacaagtacgtgcatgacgatggccgagtgctacccttgggttctcttccagccctgtgggtgcaagccgggaacaacagatgatgtatgcaggagtaaaaacaggctggtgcagacagcagagctcacaagggtg
ttcgaaatccgcaccactgatgacctcactgaggcctggctccaagaaaagttgtctttcttctgtga