

ITGA4-integrin, alpha 4 (antigen CD49D, alpha 4 subunit of VLA-4 receptor) Gene

Reference: PTXBC080190

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD49D; IA4; integrin alpha-4; 269C wild type; CD49 antigen-like family member D; VLA-4 subunit alpha; alpha 4 subunit of VLA-4 receptor; antigen CD49D, alpha-4 subunit of VLA-4 receptor; integrin alpha-IV; very late activation protein 4 receptor, alpha 4		

For research use only.

GENE DESCRIPTION	integrin, alpha 4 (antigen CD49D, alpha 4 subunit of VLA-4 receptor)
GENE ID	3676
ORIGIN SPECIES	Human
ACCESSIONS	BC080190
VECTOR	pDONR223
NCBI SYMBOL	ITGA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaatctaaagcacagcagagtgactgtacataacctttaaataatgaggttaagctgactgttcattgggttgaacccaacttcatttggatgatcaaatgatgaaatgagcctgaaacgtgcatgggtggagaaaatgaacttaacttccatgttatcaacactggcaatagtat
ggctcccaatgttagtgaggaaataatggtaccacaaattcttttagccccaactgataagctgttcaacatttggatgtccagactactactggagaatgccactttgaaaattatcaagagtggtgcatlagagcagcaaaagagtgcaatgcagaccttgaaggcatagtcagttcttctgt
ccaagactgataagaggctatgtactgcataaaagctgatccacattgtttaatttcttggtaatttgggaaaatggaagtggaagagaagccagtggtcatatccaactggaagccggccatccattttagaatggatgagacttcagcactcaagtttgaataagagcaacaggttt
tccagagccaaatccaagagtaattgaactaaacaggatgagaatgttgcgcatttctactggaaggactacatcaaaagacccaacgttattcaccatagtgattattcaagtagcttgcacttattgactctgttgcattcatatgttatgtgaaggctggcttctttaa
gacaatacaaatctatctacagaagaagaacagaagagacagttggagttatatcaacagtaaaagcaatgatgattaa