



Reference: PTXBC080554

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SFRS16; CLASP; SWAP2; CLK4-associating serine/arginine rich protein; Clk4 associating SR-related protein; splicing factor, arginine/serine-rich 16 (suppressor-of-white-apricot homolog, Drosophila); suppressor of white apricot homolog 2; CLK4 associating s		

GENE DESCRIPTION	splicing factor, arginine/serine-rich 16
GENE ID	11129
ORIGIN SPECIES	Human
ACCESSIONS	BC080554
VECTOR	pDONR223
NCBI SYMBOL	SFRS16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]