



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM198-transmembrane protein 198 Gene

Reference: PTXBC080573

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

TMEM198A; transmembrane protein 198

For research use only.

GENE DESCRIPTION	transmembrane protein 198
GENE ID	130612
ORIGIN SPECIES	Human
ACCESSIONS	BC080573
VECTOR	pDONR223
NCBI SYMBOL	TMEM198
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1