



Reference: PTXBC080574

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NGS17; TAPA; TPN; TPSN; NGS-17; TAP binding protein (tapasin); TAP-associated protein; TAP binding protein		

For research use only.

GENE DESCRIPTION	TAP binding protein (tapasin)
GENE ID	6892
ORIGIN SPECIES	Human
ACCESSIONS	BC080574
VECTOR	pDONR223
NCBI SYMBOL	TAPBP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr