

**MLLT10-myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila);
translocated to, 10 Gene**

Reference: PTXBC080577

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AF10; protein AF-10; ALL1-fused gene from chromosome 10 protein; myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila); translocated to, 10; type I AF10 protein; type III AF10 protein; type IV AF10 protein; myeloid/lymphoid or mixed-l		

For research use only.

GENE DESCRIPTION	myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila); translocated to, 10
GENE ID	8028
ORIGIN SPECIES	Human
ACCESSIONS	BC080577
VECTOR	pDONR223
NCBI SYMBOL	MLLT10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtctctagcgaccggcccggtgtcactggaggacgaggtctccatagtagatgaaggagatgattggagctgttgcgtttgctcagacgagagagctggccgagaaacccgctggtttattgcgacgggcacggcgtgcagcgtcgcggtgcatcaagcttgcattggttcaag
taccactggaccgtggtttgacggaatgtgaatctcaggagagacagccagagtgccagagctcgtctgttgcggcagctaaagtcagtggtgtgatctcagccactgcaactctgctccccgggttcaagcgattctcgtcctcagcctcccaaatggatgcaattcctgtt
ggttagcctcatctga