



Reference: PTXBC080579

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MG21; WD repeat-containing protein 13; WD repeat domain 13		
For research use only.			

GENE DESCRIPTION	WD repeat domain 13
GENE ID	64743
ORIGIN SPECIES	Human
ACCESSIONS	BC080579
VECTOR	pDONR223
NCBI SYMBOL	WDR13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]