

FCF1-FCF1 small subunit (SSU) processome component homolog (S. cerevisiae) Gene

Reference: PTXBC080600

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FCF1 rRNA-processing protein; FCF1 small subunit; rRNA-processing protein FCF1 homolog; Bka; C14orf111; CGI-35; UTP24		

For research use only.

GENE DESCRIPTION	FCF1 small subunit (SSU) processome component homolog (S. cerevisiae)
GENE ID	51077
ORIGIN SPECIES	Human
ACCESSIONS	BC080600
VECTOR	pDONR223
NCBI SYMBOL	FCF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggactgtctgtatgccaaagtgtatcccatgtataaccgattgtgtatggctgaaattgagaaattggggcagaagtatcgagtggtcttaaggattgccaaaggatccaagattgaacgattaccatgtacacacaaaggaacctatgcagatgactgcttagtacagagtaactcagcataagtggtacattgtggccacagttgaccgggacctgaaaagaagaatccgtaagattcctggagttcctcatcatgtacatttcaacctagatacaacattgaacggatgccagatgattatggagccctcgattctaa