



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZNF385C-zinc finger protein 385C Gene

Reference: PTXBC080613

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

zinc finger protein 385C; CTD-2132N18.2; CTD-2132N18.4

For research use only.

GENE DESCRIPTION	zinc finger protein 385C
GENE ID	201181
ORIGIN SPECIES	Human
ACCESSIONS	BC080613
VECTOR	pDONR223
NCBI SYMBOL	ZNF385C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaggtcagcagggggctccccggaggagcgggggcgcgcgtgtgtccagggagggtgccggacacaagccaaagagagtcacagggggccggggcgccgcgagggcccccgcctgcttccactgctgtctctcgaagtcacagtgtaattcagagaccacactgaa
cgcacatcagatcagcaggagggacacacaagaacccgtgcggggagacccccagccctccacccagcagcaccaagctcctgagacgacgacgcgtgctgtgagtatctcaagtttaaacctggccttcgaggaagcaactcaagaagcgtttggcagccgctctctgccgcgc
cgctcccacccgacgcacatgcctcgtctgtcgtcgtccggccctgccctcgcctcgcctccacacgcacatcaactctccgcctcgtccatctcggccacactgttttgcacccacagcagagctgcgcctcgtcagcagacactcgtcgtccctctttag