

LOC440461-similar to Rho GTPase activating protein 15 Gene

Reference: PTXBC080624

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to Rho GTPase activating protein 15
GENE ID	440461
ORIGIN SPECIES	Human
ACCESSIONS	BC080624
VECTOR	pDONR223
NCBI SYMBOL	LOC440461
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccggcgatgtggaggggatgtctacgtgggtggagcacccttcgagtatacccgcaaggactggcgccgctggccatctggcccaatgagcgctactggctgctggcgcgagcagcgagcactgggtggcagctggcgcgagcctggcgccaccctctacgtg
cccgcgagtagctggcgagctgcctgcctgggcaaatga