



Reference: PTXBC080636

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EPS8R2; epidermal growth factor receptor kinase substrate 8-like protein 2; EPS8-related protein 2; epidermal growth factor receptor pathway substrate 8-related protein 2; EPS8 like 2		

For research use only.

GENE DESCRIPTION	EPS8-like 2
GENE ID	64787
ORIGIN SPECIES	Human
ACCESSIONS	BC080636
VECTOR	pDONR223
NCBI SYMBOL	EPS8L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]