

LOC554249-hypothetical LOC554249 Gene

Reference: PTXBC084570

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

For research use only.

GENE DESCRIPTION	hypothetical LOC554249
GENE ID	554249
ORIGIN SPECIES	Human
ACCESSIONS	BC084570
VECTOR	pDONR223
NCBI SYMBOL	LOC554249
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaaagctgagactcagggcgagcaaccgggtcccagcggagcggcggcgcacgcggcggcacttcagcaccagtcgggtggccaccactgtgcggagatggctgcgacggtgcgcagatctcgatccaaactccctctgccagaatctggaccggaatccacca
ttgccgttttctgctgccgctggagagaatctctgaggtcccaggagagcctgctgcacggaagagatgcctcctcagfatggccgccccggagaggagcgattaa