



Reference: PTXBC089040

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	myristoylated alanine-rich protein kinase C substrate (MARCKS, 80K-L); 80K-L; MACS; PKCSL; PRKCSL; myristoylated alanine-rich C-kinase substrate; phosphomyristin; protein kinase C substrate, 80 kDa protein, light chain; myristoylated alanine rich protein		

GENE DESCRIPTION	myristoylated alanine-rich protein kinase C substrate
GENE ID	4082
ORIGIN SPECIES	Human
ACCESSIONS	BC089040
VECTOR	pDONR223
NCBI SYMBOL	MARCKS
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1