



Reference: PTXBC090852

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FIG1; LAAO; LAO; L-amino-acid oxidase; Fig-1 protein; IL4-induced protein 1; interleukin four induced 1; interleukin 4 induced 1		

For research use only.

GENE DESCRIPTION	interleukin 4 induced 1
GENE ID	259307
ORIGIN SPECIES	Human
ACCESSIONS	BC090852
VECTOR	pDONR223
NCBI SYMBOL	IL4I1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]