

OAZ1-ornithine decarboxylase antizyme 1 Gene

Reference: PTXBC093652

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AZI; OAZ; ornithine decarboxylase antizyme 1; ODC-Az; antizyme 1		

For research use only.

GENE DESCRIPTION	ornithine decarboxylase antizyme 1
GENE ID	4946
ORIGIN SPECIES	Human
ACCESSIONS	BC093652
VECTOR	pDONR223
NCBI SYMBOL	OAZ1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaaatctccctgcagcggatcctcaatagccactgcttcgagagagaagggaaggataaaccagcgccaccatccacgcccagccacatgccctcctaagcctgcacagccggcgccgagcagcagtgagaggtccagggtctccctccactgctgtagtaacc
cgggtccggggcctcgggtggtgctcctga