

C17orf67-chromosome 17 open reading frame 67 Gene

Reference: PTXBC093905

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C17orf67; chromosome 17 open reading frame 67		

For research use only.

GENE DESCRIPTION	chromosome 17 open reading frame 67
GENE ID	339210
ORIGIN SPECIES	Human
ACCESSIONS	BC093905
VECTOR	pDONR223
NCBI SYMBOL	C17orf67
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagacattgcctgtgctgctgtctcttaccttactgactgtcttctcagagacctccccgattttgacagagaagcaggccaacagctcctaagatctcggcgacaggatagaccaagcaaacccggattccccgatgagccaatgcgggaatacatgcaccacctgctgcctg
gagcatcgcctgaggagcagttcctggagcactggctgaacctcactgcgaacccactgtgacaggaacaggattcatcctgtgtaa