

OCLM-oculomedin Gene

Reference: PTXBC093987

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TISR; trabecular meshwork inducible stretch response protein		
For research use only.			

GENE DESCRIPTION	oculomedin
GENE ID	10896
ORIGIN SPECIES	Human
ACCESSIONS	BC093987
VECTOR	pDONR223
NCBI SYMBOL	OCLM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggatgtatccaccattattgttaaagatttaccttagcagacacattagtatactctctatttataaaatcctttataaaagtgtattatatggctatcctgtgtattcttcatattgttagtacttga