

RNASE3-ribonuclease, RNase A family, 3 (eosinophil cationic protein) Gene

Reference: PTXBC096060

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ECP; RAF1; RNS3; RNase 3; cytotoxic ribonuclease; ribonuclease 3; ribonuclease, RNase A family, 3; ribonuclease A family member 3		

For research use only.

GENE DESCRIPTION	ribonuclease, RNase A family, 3 (eosinophil cationic protein)
GENE ID	6037
ORIGIN SPECIES	Human
ACCESSIONS	BC096060
VECTOR	pDONR223
NCBI SYMBOL	RNASE3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttccaaaactgttcaacttccaaatttctctgttcttctgttggggttatgggtgtggagggtcactccatgccagacccccacagtttacgagggctcagtggtttgccatccagcacatcagttctgaacccccctcgatgcacattgcaatgcgggcaattaacaattatcgtatgc
gttgcaaaaacccaaatactttctctgtacaaacttttctaatgtatgttaattgttgtaaccaaaagtatacgtgccctcataacagaactcacaactgtcatcggagtagattccgggtgctttactccactgtgacctcataatccagggtcacagaatatttcaactgcaggatgc
agacagaccagggaaggagggtctatgtatgtgcatgtgacaacagagatccacgggatttccacgggtatctgtgtgttcagttcactggataccaccatctaa