



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GPR68-G protein-coupled receptor 68 Gene

Reference: PTXBC096071

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

AI2A6; GPR12A; OGR1; ovarian cancer G-protein coupled receptor 1; ovarian cancer G protein-coupled receptor, 1; sphingosylphosphorylcholine receptor; G protein-coupled receptor 68

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 68
GENE ID	8111
ORIGIN SPECIES	Human
ACCESSIONS	BC096071
VECTOR	pDONR223
NCBI SYMBOL	GPR68
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr