

SPRR2A-small proline-rich protein 2A Gene

Reference: PTXBC096108

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	small proline-rich protein 2A; SPR-2A; small proline rich protein 2A		

For research use only.

GENE DESCRIPTION	small proline-rich protein 2A
GENE ID	6700
ORIGIN SPECIES	Human
ACCESSIONS	BC096108
VECTOR	pDONR223
NCBI SYMBOL	SPRR2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttatcaacagcagcagtgcaagcagccctgccagccacctctgtgtgccccacgccaagtgtccagaacctgtccacccccgaagtgcctgagccctgccaccaccaagtgtccacagccctgccacctcagcagtgccagcagaataatctctgtgacaccttcc
ccacctgtccagtcgaagtatccaccgaagagcaagtaa