

HIST1H2BD-histone cluster 1, H2bd Gene

Reference: PTXBC096122

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2B.1B; H2B/b; H2BFB; HIRIP2; dJ221C16.6; histone H2B type 1-D; H2B histone family, member B; HIRA-interacting protein 2; histone 1, H2bd; histone H2B.1 B; histone H2B.b; histone cluster 1, H2bd; histone cluster 1 H2B family member d		

For research use only.

GENE DESCRIPTION	histone cluster 1, H2bd
GENE ID	3017
ORIGIN SPECIES	Human
ACCESSIONS	BC096122
VECTOR	pDONR223
NCBI SYMBOL	HIST1H2BD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctgaacctgctaagtcgctcgtccccaagaagggtccaagaagcgggtgactaaggctcagaagaaggacgggaagaagcgaagcgcagccgcaaggagagctattcagtgatgtgtacaaggtgctgaagcaggtccatcccgacaccggcatcttccaaggca
atggggatcatgaattccttcgtcaacgacatcttcgagcgcacgcagggcaggcttcccgctggcgcattacaacaagcgtcgcacatcacctccaggagatccagacggccgtgcgctgctgtctccggggagctggccaagcagccgtgtcggaggcaccaggcc
gtcaccagtagcaccagttccaagtaa