

HIST1H3H-histone cluster 1, H3h Gene

Reference: PTXBC096128

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H3/k; H3F1K; H3FK; histone H3.1; H3 histone family, member K; histone 1, H3h; histone H3/k; histone cluster 1, H3h; histone cluster 1 H3 family member h		

For research use only.

GENE DESCRIPTION	histone cluster 1, H3h
GENE ID	8357
ORIGIN SPECIES	Human
ACCESSIONS	BC096128
VECTOR	pDONR223
NCBI SYMBOL	HIST1H3H
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcgtacgaagcaaacagctcgaagtcaccggcgcaaggtccgcgaagcagctggccaccaaggcggtcgggaagcgtccggccaccggcggtgtcaagaagccccatcgctatcgccctgttacagtggctctcgcgagattcgccgtaccagaagtcca
ccgagctgctgatcagaagctgcttttcagcgtctggtgctgagatcgccgagacttcaagaccgacttgcgttccagagctccggtgatggcgtgcaagagcgatcgaggcctactggtgggctctttgaggacaccaacctgtgcgccatccacgccaagcgggtg
actatcatgccaagacatccagctcgacgtcgtatccgcggcgagagggctga