

TNP2-transition protein 2 (during histone to protamine replacement) Gene

Reference: PTXBC096135

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TP2; nuclear transition protein 2; TP-2; transition protein 2 (during histone to protamine replacement); transition protein 2		

For research use only.

GENE DESCRIPTION	transition protein 2 (during histone to protamine replacement)
GENE ID	7142
ORIGIN SPECIES	Human
ACCESSIONS	BC096135
VECTOR	pDONR223
NCBI SYMBOL	TNP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacacccagactcacagcttctatcaccacactcagctccatagcaacttcagccccaaagcgcacctgcacccgccattgccaacaccttcagccagagttgcagacagagccatcgtggcagccggagccagagctccagccagagcccgccagccaccgcaaccca
actggagccacagctcatccggccaccagagccagagtcaccaactagtcaccaccaagcgcacaaaaagactatgaactcccaccactctccatgcggccaccatcctgcactgccgctgccccagaacagaaagaacttgaaggccaagctgaaaaagaaaaaatg
gccaagaggatccagcagggtgtacaaaaccaagacgcggagctcaggatggaatccaactaa