

IL13-interleukin 13 Gene

Reference: PTXBC096139

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	P600; interleukin-13; interleukin 13		
For research use only.			

GENE DESCRIPTION	interleukin 13
GENE ID	3596
ORIGIN SPECIES	Human
ACCESSIONS	BC096139
VECTOR	pDONR223
NCBI SYMBOL	IL13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacccgctcctcaatcctcctctgtggcactgggcctcatggcgctttgtgaccacggctcattgctcctcacttgccttggcggtttgctcccccaggccctgtgctccctctacagccctcagggagctcattgaggagctgggtcaacatcacccagaaccagaaggctccgctctg
caatggcagcatggtatggagcatcaactgacagctggcatgtactgtgcagccctggaatccctgatcaacgtgtcaggtgcagtgccatcgagaagaccagaggatgctgagcggattctgcccgcacaagggtctcagctgggcagttttccagcttgcattgccgagacaccaa
aatcgagggtggccagtttgtaaggacctgctcttacaatttaagaaacttttcgcgaggacagttcaactga