

HIST1H1E-histone cluster 1, H1e Gene

Reference: PTXBC096168

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H1.4; H1E; H1F4; H1s-4; dJ221C16.5; histone H1.4; H1 histone family, member 4; histone 1, H1e; histone H1b; histone H1s-4; histone cluster 1, H1e; histone cluster 1 H1 family member e		

For research use only.

GENE DESCRIPTION	histone cluster 1, H1e
GENE ID	3008
ORIGIN SPECIES	Human
ACCESSIONS	BC096168
VECTOR	pDONR223
NCBI SYMBOL	HIST1H1E
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgagactgcgctgccgcgccgctgccgccccctgccgagaagactcccgtgaagaagaaggcccccaagctgcaggtgcggccaagcgaagcgtctgggccccgggtgccgagctcattactaaagctgttgcgcctccaaggagcgcagcggcgatctttg
gccgctctcaagaagcgctggcagccgctggctatgacgtggagaagaacaacagccgcatcaagctgggtctcaagagcctggtgagcaaggccaccctgtgagaccaaggccaccggcgctcgggttcctcaaaactcaacaagaagcgccctctggggaagccaagc
ctaaggctaaaaaggcagcgccgccaaggccaagccaagcagcagggcgaagaagccaagggcggggcgccaccccccaagaagcgcgaagaagacccaagaagggcgaagaagccggctgcagctgctggagccaaaaagcgaagccccga
aaaaggcgaagcagccaagccaagaagcgccaagagccagcgaaggccaagcagttaaaccaaggcggttaaccaaaagaccgccaagccaagccaagccaagaagggcgagccaagaagaagtag