

MBL2-mannose-binding lectin (protein C) 2, soluble (opsonic defect) Gene

Reference: PTXBC096179

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COLEC1; HSMBPC; MBL; MBL2D; MBP; MBP-C; MBP1; MBPD; mannose-binding protein C; collectin-1; mannan-binding lectin; mannose-binding lectin (protein C) 2, soluble (opsonic defect); mannose-binding lectin 2, soluble (opsonic defect); mannose binding lectin 2		

For research use only.

GENE DESCRIPTION	mannose-binding lectin (protein C) 2, soluble (opsonic defect)
GENE ID	4153
ORIGIN SPECIES	Human
ACCESSIONS	BC096179
VECTOR	pDONR223
NCBI SYMBOL	MBL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccctgtttccatcactccctctctctctgagatggtggcagcgtcttactcagaaactgtgacctgtgagatgcccaaaagacctgccctgcagtgtgctgtgagctctccaggcatcaacggcttcccaggcaaaagatggcggtgatggcaccaggagagaaaggggga
accaggccaagggtcagaggcttacagggtcccccctggaaagtggggtccaggaaatccagggtctctgggtcaccaggaccaaaagggtcccaaaaggagacctggaaaaagtcggatggtgatagctggctgctcagaaagaaagctctgcaaacagaaatg
gcaggtatcaaaagtggctcactctctctgggcaacaagtgggaacaagttctcctgaccaatggtgaataatgacctttgaaaaagtgaaggccttgtgtgcaagttccaggcctctgtggccacccccagggaatgctgcagagaatggagccattcagaatctcatcaaggag
gaagccttctgggcatcactgatgagaagacagaaggcagttgtgatctgacaggaaatagactgacctacacaaactggaacgaggtgaaccaacaatgctggttctgatgaagattgtgtattgctactgaaaatggccagtggaatgacgtcccctgctccacctccatct
ggccgtctgtgagttccctatctga