

PRTN3-proteinase 3 Gene

Reference: PTXBC096183

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACPA; AGP7; C-ANCA; CANCA; MBN; MBT; NP-4; NP4; P29; PR-3; PR3; myeloblastin; C-ANCA antigen; Wegener granulomatosis autoantigen; azurophil granule protein 7; leukocyte proteinase 3; neutrophil proteinase 4; serine proteinase, neutrophil; wegner autoanti		

For research use only.

GENE DESCRIPTION	proteinase 3
GENE ID	5657
ORIGIN SPECIES	Human
ACCESSIONS	BC096183
VECTOR	pDONR223
NCBI SYMBOL	PRTN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcaccggccccagccctgcccctggcgtccgtgctgctggccttgcctgagcgggtgctgcccagctgctggagatcgtggcggggcacgaggcgcagccacactcccggccctacatggcctcctgcagatgcgggggaacccgggcagccacttctcggaggca
ccttgatcccccagcttcgtgctgacggccgcgcactgctgcccacatccccagcgcctgggtgaacgtggtgctcggagcccacaacgtgctggacgagggagcccccagcagcacttctcgggtgctcaggtgtttctgaacaactacagcgggagaacaaactgaacg
acgttctcctcatccagctgagcagccagccaactcagtgctcctcgtgcccagctcagctgcccagcagggaccagccagtgccccacggcaccagtgccctggccatggcctggcgccgctgggtgcccagccccccagccagctcctgcaggagctcaatgtcacc
gtggtcaccttcttccggccacataacattgcacttctgctcctcggcgaaggccggcatctgctggagactcaggtgccccctgatctgtgatggcatcatcaaggaaatagactccttctgctgctgggatgtgccaccgccttttcctgacttctacgcggtagccctc
tacgtggactggatccgttcacgctgcgcctgtgaggccaaggccgccttga