

HOXB1-homeobox B1 Gene

Reference: PTXBC096192

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HCFP3; HOX2; HOX2I; Hox-2.9; homeobox protein Hox-B1; homeobox protein Hox-2I; homeobox B1		

For research use only.

GENE DESCRIPTION	homeobox B1
GENE ID	3211
ORIGIN SPECIES	Human
ACCESSIONS	BC096192
VECTOR	pDONR223
NCBI SYMBOL	HOXB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactataatagatgaactcctcttagagtaccactctgtaccgggggaccagcgctacagcgccacagcgccacagcgcccaactctcttcccccaagctcggctcaggcggttgacagctatgcaagcgaggccgctacgggtgggggctgtccagccctgcgttc
cagcagaactccggctatcccggccagcgcgcttcgacctgggggtgccccctccagctctgcgcctcggggtatgctcctgcgcctgcagccccagctacgggctctcagtaactacctctgggtcattcagaagagacggaggctatttcacctcgcagctacgggg
cccagctagggggcttgcgatggctacggagcaggtggagcgggtccggggccatctcccgagcatcccccttatgggaacgagcagaccgcaagctttgcacggcctatgctgaltcctctccgaggacaaggaaacacctgccccttcagaacctaacacccccacggc
ccggaccttcgactggatgaaggttaagagaacccacccaagacaggttagggctcagatgtggccacccctctccggggcccaagcatctctgttccctgcagtgacatgtcctgggtctgggctgggttctcttcttcagcggctctgtagacacaggttaa