



Reference: PTXBC096239

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HH15; HS6ST; heparan-sulfate 6-O-sulfotransferase 1; HS6ST-1; heparan-sulfate 6-sulfotransferase; heparan sulfate 6-O-sulfotransferase 1		

For research use only.

GENE DESCRIPTION	heparan sulfate 6-O-sulfotransferase 1
GENE ID	9394
ORIGIN SPECIES	Human
ACCESSIONS	BC096239
VECTOR	pDONR223
NCBI SYMBOL	HS6ST1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]