

CDX1-caudal type homeobox 1 Gene

Reference: PTXBC096252

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	caudal-type homeobox protein CDX1; homeobox protein CDX-1; caudal type homeobox transcription factor 1; caudal type homeobox transcription factor 1; caudal-type homeobox protein 1; caudal type homeobox 1		

For research use only.

GENE DESCRIPTION	caudal type homeobox 1
GENE ID	1044
ORIGIN SPECIES	Human
ACCESSIONS	BC096252
VECTOR	pDONR223
NCBI SYMBOL	CDX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtctggacaaggattcgccgggacacacgtctcgcgggagcgcagagggcggagccctacgagtgatgcggcgagcgtggcgccggaggcgccggtggcagcggttaagactcggaaccaaggacaagtaccgctgtgtctacaccgaccaccaacgcctggagctggagaaggagtttcattacagccgtttacatcacaatccggcggaatcagagctggctgccaatctggggctcactgaacggcaggtgaagatctgtttcctcaaaacggcggaaggagcgcaaaagtgaacaagaagaacagcagcagcaacagccccacagcccgatggccacgacatcacggccacccagccggccatccctggggggcctgtgtgccagcaacaccagcctcctggccacctcctccaatgcctgtgaaagaggagtttctgcatag