



Reference: PTXBC096288

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GCMA; hGCMA; chorion-specific transcription factor GCMA; GCM motif protein 1; glial cells missing homolog 1		

For research use only.

GENE DESCRIPTION	glial cells missing homolog 1 (Drosophila)
GENE ID	8521
ORIGIN SPECIES	Human
ACCESSIONS	BC096288
VECTOR	pDONR223
NCBI SYMBOL	GCM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]