

ADIPOQ-adiponectin, C1Q and collagen domain containing Gene

Reference: PTXBC096308

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACDC; ACRP30; ADIPQTL1; ADPN; APM-1; APM1; GBP28; 30 kDa adipocyte complement-related protein; adipocyte complement-related 30 kDa protein; adipose most abundant gene transcript 1 protein; adipose specific collagen-like factor; gelatin-binding protein 28;		

For research use only.

GENE DESCRIPTION	adiponectin, C1Q and collagen domain containing
GENE ID	9370
ORIGIN SPECIES	Human
ACCESSIONS	BC096308
VECTOR	pDONR223
NCBI SYMBOL	ADIPOQ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgttgcggagctgttctactgtattagctctgcccgatgaccaggaaaccacgactcaaggcccgagtcctgttccctgcccaaggggcctgcacaggttgatggcgggcatcccaggccatccggccataatggggccccagggccgtgatggcagagatg
gcacccctggtagaagggtgagaaaggagatccaggcttattgtcttaaggagacatcgggtgaaaccggagtagccggggctgaaggctcccgaggcttccgggaatccaaggcaggaaaggagaaacctggagaagggtgcctatgtataccgtcagcattcagtggtggatt
ggagacttacgttactatcccaacatgccattcgtttaccaagatcttctacaatcagcaaaaccactatgatggctccactggttaattccactgcaacattcctggcgtgactactttgcctaccacatcacagictatatgaaggatgtgaaggtcagccttccaagaaggacaaggct
atgctttcacctatgatcagtagcaggaaaataatgtgaccaggccctccgctctgtctcctgcactctgaggtggcgaccaagtctgctccaggtgtatggggaaggagagcgtaatggactctatgtgataatgacaatgactccaccttcacaggctttctctacatgacac
caactga