



Reference: PTXBC096353

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SHFM6; STHAG8; WNT-12; protein Wnt-10b; WNT-10B protein; wingless-type MMTV integration site family, member 10B; Wnt family member 10B		

For research use only.

GENE DESCRIPTION	wingless-type MMTV integration site family, member 10B
GENE ID	7480
ORIGIN SPECIES	Human
ACCESSIONS	BC096353
VECTOR	pDONR223
NCBI SYMBOL	WNT10B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]