

CDC42SE2-CDC42 small effector 2 Gene

Reference: PTXBC096703

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPEC2; CDC42 small effector protein 2; non-kinase Cdc42 effector protein SPEC2; small effector of CDC42 protein 2; CDC42 small effector 2		

For research use only.

GENE DESCRIPTION	CDC42 small effector 2
GENE ID	56990
ORIGIN SPECIES	Human
ACCESSIONS	BC096703
VECTOR	pDONR223
NCBI SYMBOL	CDC42SE2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgaattctgtgtgtttcaactgctgtattgcagaacagcctcagcctaaaaggcgacggcggtgacagaagtatgattggagagcccaaaactttgtcatagcgtcatgttgatcaggagacgtgtcagtgaatgaattcagttagctccattcagaacaaatgcagtc
caaggagggttatggaggtggaatgcctgccaatgtccagatgcagctcgttgatacgaaggcgggatag