

HIST2H2AA4-histone cluster 2, H2aa4 Gene

Reference: PTXBC096705

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2A/R; histone H2A type 2-A; histone 2, H2aa4; histone H2A/r; histone cluster 2, H2aa4; histone cluster 2 H2A family member a4		

For research use only.

GENE DESCRIPTION	histone cluster 2, H2aa4
GENE ID	723790
ORIGIN SPECIES	Human
ACCESSIONS	BC096705
VECTOR	pDONR223
NCBI SYMBOL	HIST2H2AA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggcgctgtaagcaggaggcaaggcccgcccaaggccaagtcgcgtcgtcccgctggccttcagttcccgtagggcgagtcgctgctgcgcaaggcaactacgaggagcgagtgaggcgccgcccgtctacatggctgcggtcctcgagtatct
gaccgctgagatcctggagctggcgccgaacgcggctcgggacaacaagaagacgcgcacatccctcgtcacctccagctggccatccgcaacgacgaggaaactgaacaagctgctgggcaaagtcacatcgcgccaggcgcgctcttgccaaatccaggccgctactgctcc
ctaagaagacggagagtcaccacaaggcaagggaagtga