



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RPL23AP13-ribosomal protein L23a pseudogene 13 Gene

Reference: PTXBC096706

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

For research use only.

GENE DESCRIPTION	ribosomal protein L23a pseudogene 13
GENE ID	56969
ORIGIN SPECIES	Human
ACCESSIONS	BC096706
VECTOR	pDONR223
NCBI SYMBOL	RPL23AP13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atfgcaccgcgaagtgtgaagaagggaagctctctggccgcctaaagctgaagccaaagcaaaagcgtttaaaggccaagaaggtagtgttgaagggtgtccacggccacaaaaaagaagatccgcgatgtcaccacccttcagcggcccaagacatgagactctggaggccgccga
atattctcgggaagacaccccacgggaagaacaaagctgacacctatgctatcaatgttctctcagtgagttgtccatgaagaagataaaagacacaacacccctgtgttacctgtgagtgttaaagccacaagccagcatcaaacaggctgtggaagaagctctgtgacattgat