

PPP1R14D-protein phosphatase 1, regulatory (inhibitor) subunit 14D Gene

Reference: PTXBC096716

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPI17-like; GBPI-1; GBPI1; protein phosphatase 1 regulatory subunit 14D; PKC-dependent PP1 inhibitory protein; gastrointestinal and brain-specific PP1-inhibitory protein 1; gut and brain phosphatase inhibitor 1; protein phosphatase 1 regulatory inhibitor		

For research use only.

GENE DESCRIPTION	protein phosphatase 1, regulatory (inhibitor) subunit 14D
GENE ID	54866
ORIGIN SPECIES	Human
ACCESSIONS	BC096716
VECTOR	pDONR223
NCBI SYMBOL	PPP1R14D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgtcttcaagccctgcttctgcacatctcccagccagatggggagaacccatgtaagaaggtccactgggtcttggaggagaaggacatcatccacagactcagagtcacagtcacccggactctccaagataccaggtccggagaccagccgctgacagtgaagatgaccggggccagctccagcgtggctggagatggagcaatgggtggatgctcaagttcaggagctcttccaggatcaagcaacccctctgagcctgagattgacctggaagctctcatggatctatccacagaggagcagaagactcagctggaggccattcttgggaactgcccgcgccccacagaggctttatctctgagctgctcagtcactcaagaactccggagactcagccggcctcagaataa