

OBP2B-odorant binding protein 2B Gene

Reference: PTXBC096717

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LCN14; OBPIIb; odorant-binding protein 2b; odorant-binding protein IIb; odorant binding protein 2B		

For research use only.

GENE DESCRIPTION	odorant binding protein 2B
GENE ID	29989
ORIGIN SPECIES	Human
ACCESSIONS	BC096717
VECTOR	pDONR223
NCBI SYMBOL	OBP2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaccctgttctcgtggtgtcacgtcggcctggcctgtccctgtccctcacctggaggaggatcacaggacgtggtacgtgaagccatggtgtcgataaggacttccggaggacaggaggccaggagggtgtcccagtgagggtgacagccctggcggtg
ggaagttggaagccaggttcacctcatgaggagatcggtgcatccagaagaaatcctgafgcggaagacggaggagcctggcaatacagcgccatggtggcgaggagctcatgtacgtgaggagctgccaggaggaccactacatctttactgcaagaccagcacc
atgggggcctgctccatgggaagctgtgggtaggaattctgataccaaccggaggccctggaagaatttaagaattgtgtgagcgcaaggagactctcgaggaggacattttcacgccctgcagacgggaagctgcgttccgaacactag